


Anti-Cathepsin B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07858

Item Number: CM0003780



Introduction

Anti-Cathepsin B Rabbit Monoclonal Antibody (P07858) for WB, IHC-P, ELISA. Detects human Cathepsin B, a lysosomal protease linked to tumor invasion and protein turnover.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin B Rabbit mAb
Remarks/Alias	KWE; APPS; CPSB; RECEUP; Cathepsin B
Species	Human
GeneID (Human)	1508
GeneID	1508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-156 of human Cathepsin B (P07858).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P07858
Protein Weight	38kDa
Shipping	Ice bag