


Anti-Cathepsin D Polyclonal Antibody For Wb, Ihc-P, Elisa - P07339

Item Number: CM0013394



Introduction

Rabbit polyclonal antibody against human Cathepsin D (CTSD, CPSD). Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for studies of lysosomal proteolysis, Alzheimer's disease, and breast cancer.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin D Rabbit pAb
Remarks/Alias	CPSD, CLN10, HEL-S-130P, Cathepsin D
Species	Human
Gene ID (Human)	1509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-249 of human Cathepsin D (NP_001900.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07339
Protein Weight	45kDa
Shipping	Ice bag