

Anti-Cathepsin D Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07339

Item Number: CM0002899



Introduction

Rabbit monoclonal antibody targeting human Cathepsin D (P07339), validated for WB, IHC-P, ELISA. Detects human and mouse Cathepsin D, a lysosomal aspartyl protease implicated in breast cancer and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin D Rabbit mAb
Remarks/Alias	CPSD; CLN10; HEL-S-130P; Cathepsin D
Species	Human
Gene ID (Human)	1509
Gene ID	1509
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07339
Protein Weight	45kDa
Shipping	Ice bag