


Anti-Cathepsin D Rabbit Polyclonal Antibody For Wb, Elisa - P07339

Item Number: CM0015747



Introduction

Rabbit polyclonal antibody targeting human Cathepsin D, suitable for WB and ELISA. Detects human and mouse Cathepsin D. Relevant for research in protein degradation, Alzheimer disease, and cancer.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin D Rabbit pAb
Remarks/Alias	CPSD; CLN10; HEL-S-130P; Cathepsin D
Species	Human
GeneID (human)	1509
GeneID	1509
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 65-412 of human Cathepsin D (NP_001900.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07339
Protein Weight	45kDa
Shipping	Ice bag