

Anti-Cathepsin G Rabbit Monoclonal Antibody For Wb, Elisa - P08311

Item Number: CM0003693



Introduction

Rabbit monoclonal antibody against human Cathepsin G, validated for Western Blot (WB) and ELISA. Recognizes the 29 kDa protease involved in antimicrobial activity and immune cell activation. Designed for research in immunology and infectious diseases.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin G Rabbit mAb
Remarks/Alias	CG; CATG; Cathepsin G
Species	Human
Gene ID (Human)	1511
Gene ID	1511
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 156-255 of human Cathepsin G (NP_001902.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08311
Protein Weight	29kDa
Shipping	Ice bag