

## Anti-Cathepsin G Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08311

Item Number: CM0014155



### Introduction

High-quality rabbit anti-cathepsin G polyclonal antibody validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying serine protease functions in inflammation, infection, and immunity.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cathepsin G Rabbit pAb                                                                                                   |
| Remarks/Alias    | CG; CATG; Cathepsin G                                                                                                    |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 1511                                                                                                                     |
| GeneID           | 1511                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 75-190 of human Cathepsin G (NP_001902.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IHC-P, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | P08311                                                                                                                   |
| Protein Weight   | 29kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |