

Anti-Cathepsin H Rabbit Monoclonal Antibody For Wb, Elisa - P09668

Item Number: CM0007093



Introduction

High-quality rabbit monoclonal antibody targeting human pro-cathepsin H (UniProt P09668), validated for WB and ELISA with cross-reactivity to mouse. Ideal for studying lysosomal protein degradation and cysteine protease function.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin H Rabbit mAb
Remarks/Alias	ACC4; ACC5; CPSB; ACC-4; ACC-5; Cathepsin H
Species	Human
Gene ID	1512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-250 of human Cathepsin H (NP_004381.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09668
Protein Weight	37kDa
Shipping	Ice bag