


Anti-Cathepsin K (Ctsk) Polyclonal Antibody For Wb And Elisa - P43235

Item Number: CM0020961



Introduction

Rabbit polyclonal antibody against human CTSK (Cathepsin K), validated for WB and ELISA. Reacts with human and mouse. Targets a thiol protease involved in bone resorption and matrix degradation. Ideal for osteoclast and bone remodeling research.

[Learn More](#)

Parameter	Value
Product Name	CTSK Rabbit pAb
Remarks/Alias	CTSO; PKND; PYCD; CTS02; CTS01; CTS02; CTSK
Species	Human
GeneID (Human)	1513
GeneID	1513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 115-329 of human CTSK (NP_000387.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P43235
Protein Weight	37kDa
Shipping	Ice bag