

Anti-Cathepsin L/V/K/H Polyclonal Antibody For Wb, Elisa - P07711

Item Number: CM0018559



Introduction

Rabbit polyclonal antibody against human Cathepsin L/V/K/H, validated for WB and ELISA. Reacts with human and rat. Suitable for studying lysosomal degradation, antigen processing, and coronavirus entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin L/V/K/H Rabbit pAb
Remarks/Alias	MEP; CATL; CTS1; Cathepsin L/V/K/H
Species	Human
Gene ID (Human)	1514
Gene ID	1514
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P07711
Protein Weight	38kDa
Shipping	Ice bag