

Anti-Cathepsin L/V/K/H Polyclonal Antibody For Wb, If/Icc, Elisa - P07711

Item Number: CM0016372



Introduction

Rabbit polyclonal antibody against human Cathepsin L/V/K/H, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studies on lysosomal proteolysis, viral entry, and cellular differentiation.

[Learn More](#)

Parameter	Specification
Product Name	Cathepsin L/V/K/H Rabbit pAb
Remarks/Alias	MEP; CATL; CTSL1; Cathepsin L/V/K/H
Species	Human
GeneID (Human)	1514
GeneID	1514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 114-333 of human Cathepsin L/V/K/H (NP_001903.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07711
Protein Weight	38kDa
Shipping	Ice bag