

## Anti-Cathepsin L/V/K/H Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07711/O60911/P43235/P09668

Item Number: CM0008338



### Introduction

Rabbit monoclonal antibody targeting Cathepsin L/V/K/H, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Detects a band at ~38 kDa. Lysosomal cysteine protease involved in protein degradation and antigen processing.

[Learn More](#)

| Parameter        | Value                                                                                                                              |
|------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cathepsin L/V/K/H Rabbit mAb                                                                                                       |
| Remarks/Alias    | MEP; CATL; CTS1; Cathepsin L/V/K/H                                                                                                 |
| Species          | Human                                                                                                                              |
| Gene ID (Human)  | 1514                                                                                                                               |
| Gene ID          | 1514/1515/1513/1512                                                                                                                |
| Immunogen        | Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Cathepsin L/V/K/H (P07711). |
| Source           | Rabbit                                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                  |
| SWISS            | P07711/O60911/P43235/P09668                                                                                                        |
| Protein Weight   | 38kDa                                                                                                                              |
| Shipping         | Ice bag                                                                                                                            |