


Anti-Cathepsin S Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - O70370

Item Number: CM0028525



Introduction

Anti-Cathepsin S Rabbit Polyclonal Antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA in mouse and rat. Cathepsin S is a lysosomal thiol protease involved in MHC class II antigen presentation and itch signaling. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin S Rabbit pAb
Remarks/Alias	Cats
Species	Mouse
GeneID (human)	1520
GeneID	13040
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-340 of mouse Cathepsin S (NP_067256.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O70370
Protein Weight	38 kDa
Shipping	Ice bag