


Anti-Cation-Independent M6Pr (Igf2R) Monoclonal Antibody For Wb, If/Icc, Elisa - P11717

Item Number: CM0005116



Introduction

Rabbit monoclonal antibody against human Cation-independent M6PR (IGF2R) for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for lysosomal trafficking and IGF2 signaling research.

[Learn More](#)

Parameter	Value
Product Name	Cation-independent M6PR (IGF2R) Rabbit mAb
Alias	MPR1; MPRI; CD222; CIMPR; M6P-R; MPR300; CI-M6PR; MPR 300; M6P/IGF2R; Cation-independent M6PR (IGF2R)
Species	Human
Gene ID	3482
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 2392-2491 of human Cation-independent M6PR (IGF2R) (NP_000867.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11717
Protein Weight	274 kDa
Shipping	Ice bag