


Anti-Cation-Independent M6Pr (Igf2R) Polyclonal Antibody For Wb, If/Icc, Elisa - P11717

Item Number: CM0018199



Introduction

Rabbit polyclonal antibody targeting human Cation-independent M6PR (IGF2R, CD222). Validated for WB, IF/ICC, and ELISA with reactivity across human, mouse, and rat samples. Suitable for lysosomal trafficking and IGF signaling research.

[Learn More](#)

Parameter	Value
Product Name	Cation-independent M6PR (IGF2R) Rabbit pAb
Remarks/Alias	MPR1; MPRI; CD222; CIMPR; M6P-R; MPR300; CI-M6PR; MPR 300; M6P/IGF2R; Cation-independent M6PR (IGF2R)
Species	Human
Gene ID (Human)	3482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2327-2491 of human Cation-independent M6PR (IGF2R) (NP_000867.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11717
Protein Weight	274kDa
Shipping	Ice bag