


## Anti-Cation-Independent M6Pr (Igf2R) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11717

Item Number: CM0002801



### Introduction

Rabbit monoclonal antibody targeting human Cation-independent M6PR (IGF2R). Validated for WB, IHC-P, ELISA applications with cross-reactivity for human and rat. Suitable for lysosomal trafficking and IGF2 signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                                           |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cation-independent M6PR (IGF2R) Rabbit mAb                                                                                                      |
| Remarks/Alias    | MPR1; MPRI; CD222; CIMPR; M6P-R; MPR300; CI-M6PR; MPR 300; M6P/IGF2R; Cation-independent M6PR (IGF2R)                                           |
| Species          | Human                                                                                                                                           |
| GeneID (Human)   | 3482                                                                                                                                            |
| GeneID           | 3482                                                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 2140-2300 of human Cation-independent M6PR (IGF2R) (NP_000867.3). |
| Source           | Rabbit                                                                                                                                          |
| Category         | Monoclonal Antibodies                                                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                                                |
| Cross Reactivity | Human, Rat                                                                                                                                      |
| SWISS            | P11717                                                                                                                                          |
| Protein Weight   | 274kDa                                                                                                                                          |
| Shipping         | Ice bag                                                                                                                                         |