


Anti-Cation-Independent M6Pr (Igf2R) Rabbit Polymab® For Wb, If/Icc, If-P, Ihc-P, Elisa - P11717

Item Number: CM0006334



Introduction

High-quality anti-Cation-independent M6PR (IGF2R) rabbit recombinant monoclonal antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for lysosomal trafficking and IGF2 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Cation-independent M6PR (IGF2R) Rabbit PolymAb®
Remarks/Alias	MPR1; MPRI; CD222; CIMPR; M6P-R; MPR300; CI-M6PR; MPR 300; M6P/IGF2R
Species	Human
GeneID (Human)	3482
GeneID	3482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2140-2300 of human Cation-independent M6PR (IGF2R) (NP_000867.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11717
Protein Weight	274kDa
Shipping	Ice bag