


Anti-Caveolin-1 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q03135

Item Number: CM0008957



Introduction

Mouse monoclonal antibody targeting human Caveolin-1 (VIP21, CGL3). Validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Caveolin-1 is a key structural component of caveolae membranes, involved in signal transduction and lipid transport.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-1 Mouse mAb
Remarks/Alias	CGL3; PPH3; BSCL3; LCCNS; VIP21; MSTP085
Species	Human
GeneID (human)	857
GeneID	857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Caveolin-1 (NP_001744.2).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03135
Protein Weight	17 kDa/20 kDa
Shipping	Ice bag