


Anti-Caveolin-1 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Ip, Elisa - Q03135

Item Number: CM0007032



Introduction

CamelBio offers recombinant Anti-Caveolin-1 Rabbit mAb for WB, IF/ICC, IHC-P, IP, ELISA, with cross-reactivity to Human, Mouse, and Rat. Supports caveolae formation, T-cell activation, and Wnt pathway studies.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-1 Rabbit mAb
Remarks/Alias	CGL3; PPH3; BSCL3; LCCNS; VIP21; MSTP085; Caveolin-1
Species	Human
GeneID (Human)	857
GeneID	857
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Caveolin-1 (NP_001744.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03135
Protein Weight	20kDa
Shipping	Ice bag