


## Anti-Caveolin-1 Rabbit Mab (Ko Validated) For Wb, Ihc-P, Elisa - Q03135

Item Number: CM0002869



### Introduction

High specificity anti-Caveolin-1 rabbit monoclonal antibody, KO validated, suitable for WB, IHC-P, and ELISA. Detects human, mouse, rat Caveolin-1. Ideal for caveolae biology and signaling pathway research.

[Learn More](#)

| Parameter        | Value                                                                             |
|------------------|-----------------------------------------------------------------------------------|
| Product Name     | [KO Validated] Caveolin-1 Rabbit mAb                                              |
| Remarks/Alias    | CGL3; PPH3; BSCL3; LCCNS; VIP21; MSTP085; Caveolin-1                              |
| Species          | Human                                                                             |
| Gene ID (Human)  | 857                                                                               |
| Gene ID          | 857                                                                               |
| Immunogen        | Recombinant fusion protein containing a.a. 1-80 of human Caveolin-1 (NP_001744.2) |
| Source           | Rabbit                                                                            |
| Category         | Monoclonal Antibodies                                                             |
| Application      | WB, IHC-P, ELISA                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                 |
| SWISS            | Q03135                                                                            |
| Protein Weight   | 17kDa/20kDa                                                                       |
| Shipping         | Ice bag                                                                           |