

Anti-Caveolin-2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P51636

Item Number: CM0011905

Introduction

Rabbit monoclonal antibody targeting human Caveolin-2, validated for WB, IHC-P, IP, and ELISA. Reacts with human and mouse, detecting an 18kDa caveolar scaffolding protein involved in MAPK signaling and caveolae formation.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-2 Rabbit mAb
Remarks/Alias	CAV; Caveolin-2
Species	Human
GenelD (Human)	858
GenelD	858
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51636
Protein Weight	18kDa
Shipping	Ice bag