

Anti-Caveolin-2 Rabbit Pab For Wb - P51636

Item Number: CM0021000

Introduction

High-quality anti-Caveolin-2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Recombinant immunogen. Ideal for caveolae formation and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-2 Rabbit pAb
Remarks/Alias	CAV; Caveolin-2
Species	Human
GeneID (Human)	858
GeneID	858
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human CAV2 (NP_001224.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51636
Protein Weight	18kDa
Shipping	Ice bag