

Anti-Caveolin-3 Polyclonal Antibody For Wb, Ihc-P, Elisa - P56539

Item Number: CM0030639

Introduction

Rabbit polyclonal antibody against human Caveolin-3, validated for Western blot, IHC-P, and ELISA. Cross-reacts with human and mouse (predicted band ~17kDa). For studying sarcolemma repair, caveolae, and muscle diseases.

[Learn More](#)

Parameter	Specification
Product Name	Caveolin-3 Rabbit pAb
Remarks/Alias	LQT9; MPDT; RMD2; VIP21; LGMD1C; VIP-21; Caveolin-3
Species	Human
Gene ID (Human)	859
Gene ID	859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-151 of human Caveolin-3 (NP_001225.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56539
Protein Weight	17kDa
Shipping	Ice bag