

Anti-Caveolin-3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P56539

Item Number: CM0009449

Introduction

High-quality rabbit monoclonal antibody targeting human Caveolin-3 (P56539). Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for investigating caveolar function and muscle-specific diseases.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-3 Rabbit mAb
Remarks/Alias	LQT9; MPDT; RMD2; VIP21; LGMD1C; VIP-21; Caveolin-3
Species	Human
GeneID (Human)	859
GeneID	859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-151 of human Caveolin-3 (P56539)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56539
Protein Weight	17kDa
Shipping	Ice bag