

Anti-Cavin2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95810

Item Number: CM0021132

Introduction

High-quality anti-CAVIN2 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for investigating caveolae-associated protein 2 (cavin-2/PS-p68) function in membrane curvature and caveolae biology.

[Learn More](#)

Parameter	Value
Product Name	CAVIN2 Rabbit pAb
Remarks/Alias	SDR; SDPR; PS-p68; cavin-2; CAVIN2
Species	Human
Gene ID (Human)	8436
Immunogen	Synthetic Peptide corresponding to a sequence within amino acids 325-425 of human CAVIN2 (NP_004648.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95810
Protein Weight	47KDa
Shipping	Ice bag