

Anti-Cbfa2T2 Polyclonal Antibody For Wb, If/Icc, Elisa - 043439

Item Number: CM0030037

Introduction

Rabbit polyclonal antibody against human CBFA2T2, a transcriptional corepressor involved in ESC pluripotency, for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	CBFA2T2 Rabbit pAb
Remarks/Alias	EHT; p85; MTGR1; ZMYND3
Species	Human
GeneID (Human)	9139
GeneID	9139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 325-604 of human CBFA2T2 (NP_005084.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43439
Protein Weight	67kDa
Shipping	Ice bag