


Anti-Cbfa2T3 Rabbit Polyclonal Antibody For Wb, Elisa - O75081

Item Number: CM0015334



Introduction

Rabbit polyclonal antibody against transcriptional corepressor CBFA2T3 (ETO2/MTG16) for WB and ELISA. Validated in human, mouse, and rat samples. Ideal for transcription regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CBFA2T3 Rabbit pAb
Remarks/Alias	ETO2; MTG16; MTGR2; ZMYND4; RUNX1T3; CBFA2T3
Species	Human
Gene ID (Human)	863
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 593-653 of human CBFA2T3 (NP_005178.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75081
Protein Weight	71kDa
Shipping	Ice bag