

Anti-Cbfb Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q13951

Item Number: CM0012025



Introduction

High-quality rabbit monoclonal antibody against human CBFB (Core-binding factor subunit beta). Validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying CBF/RUNX transcription factor complexes.

[Learn More](#)

Parameter	Value
Product Name	CBFB Rabbit mAb
Remarks/Alias	CLCD2; PEBP2B
Species	Human
Gene ID (Human)	865
Gene ID	865
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human CBFB (NP_001746.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13951
Protein Weight	22 kDa
Shipping	Ice bag