


Anti-Cbl Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P22681

Item Number: CM0013972



Introduction

High-quality rabbit polyclonal antibody against human CBL (E3 ubiquitin-protein ligase) for WB, IHC-P, IF/ICC, IP, and ELISA. Validated in human, mouse, and rat tissues. Ideal for studying ubiquitination, signaling, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	CBL Rabbit pAb
Remarks/Alias	CBL2; NSLL; C-CBL; RNF55; FRA11B; CBL
Species	Human
Gene ID (Human)	867
Gene ID	867
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 850-906 of human CBL (NP_005179.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22681
Protein Weight	100kDa
Shipping	Ice bag