


Anti-Cbl Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P22681

Item Number: CM0007514



Introduction

High-quality rabbit monoclonal antibody targeting human E3 ubiquitin-protein ligase CBL (P22681) for Western blot, IHC-P, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for signal transduction and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CBL Rabbit mAb
Remarks/Alias	CBL2; NSLL; C-CBL; RNF55; FRA11B; CBL
Species	Human
Gene ID (Human)	867
Gene ID	867
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 850-906 of human CBL (NP_005179.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22681
Protein Weight	100kDa
Shipping	Ice bag