


Anti-Cbl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22681

Item Number: CM0025759



Introduction

CBL Rabbit Polyclonal Antibody validated in WB, IHC-P, and ELISA. Detects human CBL, cross-reacts with mouse and rat. Designed for research on E3 ubiquitin ligase signaling and receptor tyrosine kinase regulation.

[Learn More](#)

Parameter	Value
Product Name	CBL Rabbit pAb
Remarks/Alias	CBL2; NSLL; C-CBL; RNF55; FRA11B; CBL
Species	Human
GeneID (human)	867
GeneID	867
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 850-906 of human CBL (NP_005179.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22681
Protein Weight	100kDa
Shipping	Ice bag