

Anti-Cblb Rabbit Polyclonal Antibody For Wb, Elisa - Q13191

Item Number: CM0024299



Introduction

Rabbit polyclonal antibody against CBLB (Cbl-b, RNF56) validated for Western blot (WB) and ELISA. Detects human, mouse, and rat CBLB. E3 ubiquitin ligase regulating immune receptor signaling; ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CBLB Rabbit pAb
Remarks/Alias	Cbl-b; RNF56; Nbla00127; CBLB
Species	Human
GeneID (human)	868
GeneID	868
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 883-982 of human CBLB (NP_001308717.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13191
Protein Weight	109kDa
Shipping	Ice bag