


Anti-Cblc Rabbit Polyclonal Antibody For Wb, Elisa - Q9ULv8

Item Number: CM0025575



Introduction

Rabbit polyclonal antibody targeting human CBLC (E3 ubiquitin-protein ligase CBL-C), validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for ubiquitination and EGFR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CBLC Rabbit pAb
Remarks/Alias	CBL-3; RNF57; CBL-SL; CBLC
Species	Human
Gene ID (Human)	23624
Gene ID	23624
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULV8
Protein Weight	52kDa
Shipping	Ice bag