

Anti-Cbll1 Rabbit Polyclonal Antibody For Wb And Elisa - Q75N03

Item Number: CM0017644

Introduction

Rabbit polyclonal antibody targeting human CBLL1 (HAKAI/RNF188), an E3 ubiquitin ligase. Validated for Western blot and ELISA, cross-reacts with mouse. Suitable for ubiquitination and m6A methylation pathway research.

[Learn More](#)

Parameter	Value
Product Name	CBLL1 Rabbit pAb
Remarks/Alias	HAKAI; RNF188; CBLL1
Species	Human
Gene ID (Human)	79872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human CBLL1 (NP_079090.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q75N03
Protein Weight	55kDa
Shipping	Ice bag