

Anti-Cbp80 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q09161

Item Number: CM0017481



Introduction

High-quality rabbit polyclonal antibody against CBP80 (NCBP1) for Western blot, IHC-P, and ELISA. Validated in human and mouse. CBP80 is a component of the cap-binding complex involved in mRNA processing and export.

[Learn More](#)

Parameter	Value
Product Name	CBP80 Rabbit pAb
Remarks/Alias	NCBP; Sto1; CBP80
Species	Human
Gene ID	4686 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CBP80 (NP_002477.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q09161
Protein Weight	92kDa
Shipping	Ice bag