


Anti-Cbr1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P16152

Item Number: CM0007140



Introduction

CBR1 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets carbonyl reductase 1, involved in xenobiotic metabolism and prostaglandin synthesis. Ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	CBR1 Rabbit mAb
Remarks/Alias	CBR; hCBR1; PG-9-KR; SDR21C1; CBR1
Species	Human
Gene ID	873
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16152
Protein Weight	30kDa
Shipping	Ice bag