

Anti-Cbr1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P16152

Item Number: CM0002093

Introduction

Rabbit monoclonal antibody to human Carbonyl reductase 1 (CBR1). Validated for WB, IHC-P, ELISA. Reacts with human and mouse. Suitable for studying carbonyl reduction, glucocorticoid metabolism, and drug detoxification.

[Learn More](#)

Parameter	Value
Product Name	CBR1 Rabbit mAb
Remarks/Alias	CBR; hCBR1; PG-9-KR; SDR21C1; CBR1
Species	Human
GeneID (human)	873
GeneID	873
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16152
Protein Weight	30kDa
Shipping	Ice bag