


Anti-Cbr1 Rabbit Polyclonal Antibody For Wb, Elisa - P16152

Item Number: CM0027810



Introduction

Rabbit polyclonal antibody against human CBR1 (Carbonyl reductase [NADPH] 1), validated for WB and ELISA. Reacts with human and mouse samples. Useful for studying carbonyl metabolism and drug detoxification.

[Learn More](#)

Parameter	Specification
Product Name	CBR1 Rabbit pAb
Remarks/Alias	CBR; hCBR1; PG-9-KR; SDR21C1; CBR1
Species	Human
GeneID (human)	873
GeneID	873
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human CBR1 (NP_001748.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16152
Protein Weight	30kDa
Shipping	Ice bag