


Anti-Cbr1 Rabbit Polyclonal Antibody For Wb - P16152

Item Number: CM0022013



Introduction

Rabbit polyclonal antibody against human CBR1, knockout validated for Western blot and ELISA. Detects endogenous CBR1 in human and mouse samples. Investigate carbonyl reduction, drug metabolism, and prostaglandin synthesis. Molecular weight ~30 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CBR1 Rabbit pAb
Remarks/Alias	CBR; hCBR1; PG-9-KR; SDR21C1; R1
Species	Human
GeneID (human)	873
GeneID	873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human CBR1 (NP_001748.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16152
Protein Weight	30kDa
Shipping	Ice bag