

Anti-Cbr4 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N4T8

Item Number: CM0030983



Introduction

Rabbit polyclonal antibody targeting human CBR4 (SDR45C1), validated for WB and ELISA. Cross-reacts with rat. CBR4 is a mitochondrial 3-oxoacyl-[ACP] reductase involved in fatty acid biosynthesis and quinone detoxification. Ideal for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CBR4 Rabbit pAb
Remarks/Alias	SDR45C1; CBR4
Species	Human
Gene ID (Human)	84869
Gene ID	84869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human CBR4 (NP_116172.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8N4T8
Protein Weight	25kDa
Shipping	Ice bag