

Anti-Cbs Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P35520

Item Number: CM0013695



Introduction

CBS Rabbit pAb for WB, IF-P, ELISA. Detects 61 kDa cystathionine beta-synthase in human, mouse, rat. Ideal for homocysteine metabolism and hydrogen sulfide signaling research.

[Learn More](#)

Parameter	Value
Product Name	CBS Rabbit pAb
Remarks/Alias	CBSL; HIP4; CBS
Species	Human
GeneID (Human)	875
GeneID	875
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35520
Protein Weight	61kDa
Shipping	Ice bag