


## Anti-Cbx1/Hp1 Beta Monoclonal Antibody For Wb, Ihc-P, Elisa - P83916

Item Number: CM0006281



### Introduction

Rabbit monoclonal antibody targeting human CBX1/HP1 beta (HP1 $\beta$ ), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Detects a 21kDa heterochromatin protein involved in epigenetic repression. Ideal for chromatin and epigenetics research.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CBX1/HP1 beta Rabbit mAb                                                                                                 |
| Remarks/Alias    | CBX; M31; MOD1; p25beta; HP1-BETA; HP1Hsbeta; HP1Hs-beta; CBX1/HP1 beta                                                  |
| Species          | Human                                                                                                                    |
| Gene ID (Human)  | 10951                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CBX1/HP1 beta (NP_006798.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, IHC-P, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | P83916                                                                                                                   |
| Protein Weight   | 21kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |