


Anti-Cbx7 Polyclonal Antibody For Wb - 095931

Item Number: CM0022903



Introduction

Rabbit polyclonal antibody targeting human CBX7, validated for WB applications. Detects a band at 28kDa. Ideal for studying Polycomb group-mediated transcriptional repression and cancer epigenetics. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CBX7 Rabbit pAb
Remarks/Alias	CBX7
Species	Human
GeneID (Human)	23492
GeneID	23492
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 162-251 of human CBX7 (NP_783640.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	095931
Protein Weight	28kDa
Shipping	Ice bag