


Anti-Cbx8 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Hc52

Item Number: CM0030228



Introduction

Rabbit polyclonal antibody targeting CBX8 (Chromobox protein homolog 8) for WB, IHC-P, ELISA applications. Validated in human, mouse, rat. Recombinant immunogen. Ideal for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	CBX8 Rabbit pAb
Remarks/Alias	PC3; RC1; CBX8
Species	Human
Gene ID (Human)	57332
Gene ID	57332
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human CBX8 (NP_065700.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC52
Protein Weight	43kDa
Shipping	Ice bag