

Anti-Ccar1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q8Ix12

Item Number: CM0015914

Introduction

CCAR1 rabbit polyclonal antibody validated for WB, IP, ELISA in human, mouse, rat. CCAR1 is a transcriptional coregulator and apoptosis regulator involved in cell cycle progression and cancer.

[Learn More](#)

Parameter	Value
Product Name	CCAR1 Rabbit pAb
Remarks/Alias	CCAR1
Species	Human
GeneID (Human)	55749
GeneID	55749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CCAR1 (NP_060707.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IX12
Protein Weight	133kDa
Shipping	Ice bag