

Anti-Ccar2 Rabbit Polyclonal Antibody For Wb, Ihc-P, And Elisa - Q8N163

Item Number: CM0026957



Introduction

CCAR2 Rabbit Polyclonal Antibody reacts with human, mouse, and rat CCAR2 (DBC1). Validated for Western blot, IHC-P, and ELISA. Recombinant immunogen. Suitable for apoptosis, transcriptional regulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CCAR2 Rabbit pAb
Remarks/Alias	DBC1; DBC-1; NET35; p30DBC; p30 DBC; KIAA1967; CCAR2
Species	Human
Gene ID (Human)	57805
Gene ID	57805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 654-923 of human CCAR2 (NP_066997.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N163
Protein Weight	103kDa
Shipping	Ice bag