

Anti-Ccbl1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q16773

Item Number: CM0029080



Introduction

Rabbit polyclonal antibody against CCBL1 (Kynurenine--oxoglutarate transaminase 1), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Useful for tryptophan catabolism and cysteine conjugate beta-lyase studies.

[Learn More](#)

Parameter	Value
Product Name	CCBL1 Rabbit pAb
Remarks/Alias	GTK; KAT1; KATI; CCBL1
Species	Human
GeneID (human)	883
GeneID	883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CCBL1 (NP_001116144.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16773
Protein Weight	48kDa
Shipping	Ice bag