


Anti-Ccdc124 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Ct7

Item Number: CM0017235



Introduction

Rabbit polyclonal antibody targeting human CCDC124, a ribosome hibernation and cytokinesis regulator. Validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for cell cycle and translational research.

[Learn More](#)

Parameter	Value
Product Name	CCDC124 Rabbit pAb
Remarks/Alias	Lso2; oxs1; CCDC124
Species	Human
Gene ID (Human)	115098
Gene ID	115098
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human CCDC124 (NP_612451.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96CT7
Protein Weight	26kDa
Shipping	Ice bag