

Anti-Ccdc137 Polyclonal Antibody For Wb, Elisa - Q6Pk04

Item Number: CM0023106

Introduction

CCDC137 Rabbit pAb is a polyclonal antibody targeting human coiled-coil domain-containing protein 137, validated for Western blot and ELISA. Suitable for studying chromosome-associated protein interactions.

[Learn More](#)

Parameter	Value
Product Name	CCDC137 Rabbit pAb
Remarks/Alias	RaRF, CCDC137
Species	Human
GeneID (Human)	339230
GeneID	339230
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 213-288 of human CCDC137 (NP_954981.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6PK04
Protein Weight	33kDa
Shipping	Ice bag