


## Anti-Ccdc47 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96A33

Item Number: CM0014503



### Introduction

Rabbit polyclonal antibody against CCDC47, a key subunit of the PAT/MPT complex. Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying ER protein insertion and calcium homeostasis.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | CCDC47 Rabbit pAb                                                                                                   |
| Remarks/Alias    | THNS; GK001; MSTP041; CCDC47                                                                                        |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 57003                                                                                                               |
| GeneID           | 57003                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 224-483 of human CCDC47 (NP_064583.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q96A33                                                                                                              |
| Protein Weight   | 56kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |