

Anti-Ccdc50 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ivm0

Item Number: CM0018992



Introduction

CCDC50 Rabbit pAb from CamelBio is a polyclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat CCDC50 (YMER). Ideal for studies on EGFR signaling and cytoplasmic protein interactions.

[Learn More](#)

Parameter	Value
Product Name	CCDC50 Rabbit pAb
Remarks/Alias	YMER, C3orf6, DFNA44, CCDC50
Species	Human
GeneID (Human)	152137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-147 of human CCDC50 (NP_777568.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IVM0
Protein Weight	36kDa
Shipping	Ice bag